



Reliable aerosol delivery that you can trust

InnoSpire Elegance compressor nebulizer system

A powerful compressor and a compact design combine to create a nebulizer system you and your patients will value. Elegance, part of the InnoSpire family of compressor nebulizer systems from Philips Respironics, provides fast, efficient and reliable drug delivery to patients at home.

This elegant performer weighs only 3.3 pounds and works efficiently with proven SideStream aerosol technology. It includes a handset, docking station, filter and power cord and is covered by a 5-year warranty. Elegance is also compliant with the international safety and essential performance standard IEC 60601-1 third edition and its amendments.

PHILIPS
RESPIRONICS
sense and simplicity

InnoSpire Elegance



InnoSpire Elegance and SideStream technology

Elegance and SideStream disposable and reusable nebulizers combine to provide fast and efficient drug delivery. SideStream nebulizers feature an active venturi system that provides air flow in addition to that supplied by the compressor, resulting in faster drug delivery and shorter treatment times. They are easy to use with fewer parts to assemble or clean than other conventional nebulizers.



Venturi system

provides a faster delivery rate resulting in shorter treatment times than other conventional nebulizers

The breath-enhanced SideStream Plus nebulizer features an easy-action inspiratory valve that opens on inspiration to boost medication delivery and closes on exhalation, preserving medication and reducing waste.¹

Less pressure is required to open the inspiratory valve compared to other breath-enhanced nebulizers.* And there's no need to remove the valve when using a mask, making treatments easy for pediatric and elderly patients.



SideStream Plus

¹ Performance characterization of the new SideStream Plus breath enhanced jet nebulizer.

* Data on file

| Product specifications | |
|-------------------------|--------------------|
| Weight | 3.3 lbs. |
| Size | 4.2" x 6.5" x 6.5" |
| Electrical rating/power | 120V/60Hz, 1.7 Amp |
| Nebulization capacity | 8 ml |
| MMAD | (80%) <5 microns |
| Treatment time | 5-7 minutes |
| Maximum pressure | 44 psi |
| Maximum flow rate | 9.4 LPM |
| Operating flow rate | 7 LPM @ 10 psi |
| Mode of operation | Continuous |
| Warranty | 5 years |

| Ordering information | | |
|---|----------------|---------------|
| Product description | Reorder number | NDC number |
| InnoSpire Elegance with SideStream disposable nebulizer | 1100134 | 08373-0134-00 |
| InnoSpire Elegance with SideStream disposable and reusable nebulizers | 1099969 | 08373-9969-00 |
| InnoSpire Elegance with SideStream Plus | 1102828 | 08373-2828-00 |
| Replacement filter (4-pack) | 1102088 | 08373-2088-00 |
| Standard carrying case | 1101452 | 08373-0452-00 |
| Premium carrying case | 1026568 | 08373-0568-00 |
| InnoSpire Elegance compressor with SideStream reusable nebulizer DVD | 1101988 | N/A |
| InnoSpire Elegance compressor with SideStream Plus nebulizer DVD | 1101989 | N/A |

Economical

| Product description | HCPCS |
|--|-------|
| InnoSpire Elegance compressor nebulizer system | E0570 |
| SideStream disposable nebulizer | A7003 |
| SideStream reusable nebulizer | A7005 |
| SideStream Plus reusable nebulizer | A7005 |

Please visit www.philips.us/elegance for more information



Respironics, InnoSpire, Elegance, and SideStream are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts this device to sale by or on the order of a physician.

DD MCI 4104969 PN 1101728

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International
Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282

Respironics New Jersey, Inc.
5 Wood Hollow Road
Parsippany, NJ 07054
+1 800 345 6443
+1 724 387 4000
www.philips.com/respironics